



Clip Applikatoren
Clip Applicators





Clip-Applikatoren *Clip-Applicators*


①




Clip-Applikator
drehbar und zerlegbar,
Ø 5.0 mm, AL 330 mm
Clip-applicator
rotatable and detachable,
Ø 5.0 mm, WL 330 mm



②



Clip-Applikator
drehbar und zerlegbar,
Ø 10.0 mm, AL 330 mm
Clip-applicator
rotatable and detachable,
Ø 10.0 mm, WL 330 mm



③



Clip-Applikator
drehbar, Spüleinsatz, nicht zerlegbar,
Ø 10.0 mm, AL 330 mm
Clip-applicator
rotatable, cleaning port, not detachable,
Ø 10.0 mm, WL 330 mm





Clip-Applikatoren Clip-Applicators




①

Clip-Applikator
drehbar und zerlegbar,
Ø 5.0 mm, AL 330 mm
*Clip-applicator
rotatable and detachable,
Ø 5.0 mm, WL 330 mm*



②

Clip-Applikator
drehbar und zerlegbar,
Ø 10.0 mm, AL 330 mm
*Clip-applicator
rotatable and detachable,
Ø 10.0 mm, WL 330 mm*



③

Clip-Applikator
drehbar, Spüleinsatz,
nicht zerlegbar,
Ø 10.0 mm, AL 330 mm
*Clip-applicator
rotatable, cleaning port,
not detachable,
Ø 10.0 mm, WL 330 mm*

Für / for

Ethicon-Titanclips

small	17.190.01	17.202.01	17.200.01
medium	17.190.02	17.202.02	17.200.02
medium-large		17.202.03	17.200.03
large		17.202.04	17.200.04

Pilling-Titanclips

small	17.190.10	17.202.10	17.200.10
medium	17.190.11	17.202.11	17.200.11
medium-large		17.202.12	17.200.12

Weck-Titanclips

small	17.190.20	17.202.20	17.200.20
medium	17.190.21	17.202.21	17.200.21
medium-large		17.202.22	17.200.22

Horizon-Titanclips

small	17.190.30	17.202.30	17.200.30
medium	17.190.31	17.202.31	17.200.31
medium-large		17.202.32	17.200.32

Für Clips weiterer Marken und Größen können Clip-Applikatoren auf Anfrage gefertigt werden.
Upon request clip-applicators can be produced in different sizes and types.