


## Anterior Cruciate Ligament Reconstruction **PRESS FIT** Technique



## **Biological** ACL Reconstruction Operation Implant – Free / **with Bone Material**



## The **ALL PRESS FIT** technique realises the 10 rules for ideal anatomical ACL-Reconstruction:

- Fixation near the original insertion
- Broad proximal insertion for imitation the bundles
- No foreign material
- Biological healing-support by pressfit fixation
- MRI-Control without artefacts
- No bonedefects for easy revision in the same technique
- Primary stability
- Fast rehabilitation
- Secure and easy to reproduce operation-technique
- With low costs

## Main Indications

- Bone biopsies
- Harvesting of autologous bone material
- Osteocartilage Defects (Single-and composite grafts) e.g.
  - Osteochondrosis dissecans
  - Traumatil cartilageulcus
  - Unilateral Gonarthrosis
- Treatment of osteochondral defects with bone grafting in combination with matrix-induced ACT
- Osteosynthetical Treatments (Fracture Dowelling)
  - Flake fracture/Osteochondral fractures
  - Fractures of the joint with bone-cartilage-fragments for refixation
  - Fractures of the condyles
  - Head of the Femur Necroses (Avascular Necroses of the Bone)
  - Metaphyseal Fractures
  - Spondylodeses
  - Arthrodeses of Symphysis, upper and lower jump joint, Root of the foot
  - Fractures and Nonunions of the Scaphoid and Talus

| Art.Nr. /part no. | Artikelbild/picture   | Artikel-Bezeichnung/part description                                       |
|-------------------|---|--|
| 23-0100-00        |    | Tibiales Zielgerät /Tibial Target Device                                   |
| 23-0101-00        |    | Zielarm für Tibiales Zielgerät /<br>Aiming device for Tibial Target Device |
| 23-0110-00        |    | Femorales Zielgerät/Femoral Target Device Ø8-7                             |
| 23-0111-00        |  | Femorales Zielgerät/Femoral Target Device Ø9-8                             |
| 23-0112-00        |  | Femorales Zielgerät/Femoral Target Device Ø10-9                            |
| 23-0113-00        |  | Zielgerät für Doppel-K-Draht/Tibial Aimer for Double K-Wire                |
| 23-0120-00        |  | Kronenhohlfräse/Crown Cutter Ø7-6/165 mm                                   |
| 23-0121-00        |  | Kronenhohlfräse/Crown Cutter Ø8-7/165 mm                                   |