

Promisemed

VeriEcto® Co-axial Biopsy Devices

Improve Accuracy
and Efficiency

Promisemed


FEATURES AND BENEFITS

Promisemed VeriEcto® Co-axial Biopsy Devices are compatible with VeriEcto® biopsy needles, helping with accurate location and allowing multiple sampling from one puncture with minimal trauma to organs and/or surrounding tissue.



Accurate Location

- Stiff cannula and extra sharp trocar tip for straighter penetration



- Adjustable depth stopper prevents forward movement



- Centimeter depth markings facilitate accurate depth control




Multiple Sampling

- Provide a clear path that makes multiple sampling from one puncture possible to minimize trauma to organs and/or surrounding tissue and improve efficiency
- Available with trocar tip stylet and blunt tip stylet




Compatible

- Compatible with VeriEcto® biopsy needles. One gauge size larger and 30mm length shorter than corresponding Verisafe® biopsy needles with the same color



Cross-section of Needle



Needle Length

ORDERING INFORMATION

VeriEcto®	REF		Gauge Size and Needle Length	
	Co-axial Biopsy Devices	Co-axial Biopsy Devices with Blunt Tip Stylet	Co-axial Biopsy Devices (with Blunt Tip Stylet)	Compatible Biopsy Needles
11G	CBN-1107	CBN-1107R	3.0 (11G) ×70mm	2.7(12G) × 100mm
	CBN-1113	CBN-1113R	3.0 (11G) ×130mm	2.7(12G) × 160mm
13G	CBN-1307	CBN-1307R	2.4(13G) ×70mm	2.1(14G) × 100mm
	CBN-1313	CBN-1313R	2.4(13G) ×130mm	2.1(14G) × 160mm
15G	CBN-1507	CBN-1507R	1.8(15G) ×70mm	1.6(16G) × 100mm
	CBN-1513	CBN-1513R	1.8(15G) ×130mm	1.6(16G) × 160mm
17G	CBN-1707	CBN-1707R	1.4(17G) ×70mm	1.2(18G) × 100mm
	CBN-1713	CBN-1713R	1.4(17G) ×130mm	1.2(18G) × 160mm
19G	CBN-1907	CBN-1907R	1.1(19G) ×70mm	0.9(20G) × 100mm
	CBN-1913	CBN-1913R	1.1(19G) ×130mm	0.9(20G) × 160mm

※ Co-axial Biopsy Device Only. Compatible biopsy needles are not included.
※ If you have further questions or need other specifications, please contact us.



PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com