

Promisemed

VeriEcto® Fine Biopsy Needles

HARVEST CELLS
FOR
CYTOLOGICAL EXAMINATIONS



Promisemed

FEATURES AND BENEFITS

Designed for harvesting cells for cytological examinations, Promisemed VeriEcto® Fine Biopsy Needles are available with Chiba tips.



- **Extra sharp Chiba tip** to ease penetration



- **Smooth inner wall** of the needle tube to reduce trauma and harvest intact cells



- **Adjustable depth stopper** prevents forward movement and aids in locating the needle tip at the biopsy site




- **Centimeter depth markings** facilitate accurate depth control






- **Transparent needle base** for good visibility during aspiration procedure



- **Standard luer connector** for easy drug injection and apocenos



- **Color code** identifies needle gauges



ORDERING INFORMATION



Gauge Size		Needle Length	REF
16G	1.6mm	150mm	FBN-1615C
		200mm	FBN-1620C
18G	1.2mm	150mm	FBN-1815C
		200mm	FBN-1820C
19G	1.1mm	150mm	FBN-1915C
		200mm	FBN-1920C
20G	0.9mm	150mm	FBN-2015C
		200mm	FBN-2020C
21G	0.8mm	150mm	FBN-2115C
		200mm	FBN-2120C
22G	0.7mm	150mm	FBN-2215C
		200mm	FBN-2220C
23G	0.6mm	150mm	FBN-2315C
		200mm	FBN-2320C

※ If you have further questions or need other specifications, please contact us.



PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com



6666-FBN-VeriEcto-H-00 (EN)