

Promisemed

VeriEcto® Semi-automatic Biopsy Needles

Get
Ideal Samples
Meet Multiple Demands




Promisemed

FEATURES AND BENEFITS


Promisemed VeriEcto® Semi-automatic Biopsy Needles offer many of the same features that our automatic biopsy needles possess. The difference is that you need to advance the stylet manually.




Get Ideal Samples



- Smaller and quieter vibration when fired
- Echogenic tip enhances visualization under ultrasound





- Extra sharp trocar tip to ease penetration




- Extra sharp cutting cannula to minimize trauma and get ideal samples


Meet Multiple Demands

- Applicable for most organs such as breast, kidney, lung, liver, lymph gland and prostate






• Two specimen notches for your choice



Designed for small tumors and areas with rich blood vessels.



Designed for other soft tissues

• Optional VeriEcto® Co-axial Biopsy Devices enhance efficiency and accuracy



User-friendly

- Smooth stylet advancement
- Ergonomic plunger and finger grips, as well as lightweight design for comfortable and precise control



Plunger



Finger grips

• Safety button to avoid accidental triggering



ORDERING INFORMATION


Promisemed VeriEcto® Semi-automatic Biopsy Needles	
REF	 Gauge Size and Needle Length
SBN-1410	2.1 (14G) ×100mm
SBN-1416	2.1 (14G) ×160mm
SBN-1420	2.1 (14G) ×200mm
SBN-1610	1.6 (16G) ×100mm
SBN-1616	1.6 (16G) ×160mm
SBN-1620	1.6 (16G) ×200mm
SBN-1810	1.2 (18G) ×100mm
SBN-1816	1.2 (18G) ×160mm
SBN-1820	1.2 (18G) ×200mm
SBN-2010	0.9 (20G) ×100mm
SBN-2016	0.9 (20G) ×160mm
SBN-2020	0.9 (20G) ×200mm

※ If you have further questions or need other specifications, please contact us.

Promisemed VeriEcto® Semi-automatic Biopsy Needles with Co-axial Biopsy Devices		
REF	Gauge Size and Needle Length	
	 Semi-automatic Biopsy Needle	 Co-axial Biopsy Device
SBN-1410C	2.1 (14G) ×100mm	2.4 (13G) ×70mm
SBN-1416C	2.1 (14G) ×160mm	2.4 (13G) ×130mm
SBN-1420C	2.1 (14G) ×200mm	2.4 (13G) ×170mm
SBN-1610C	1.6 (16G) ×100mm	1.8 (15G) ×70mm
SBN-1616C	1.6 (16G) ×160mm	1.8 (15G) ×130mm
SBN-1620C	1.6 (16G) ×200mm	1.8 (15G) ×170mm
SBN-1810C	1.2 (18G) ×100mm	1.4 (17G) ×70mm
SBN-1816C	1.2 (18G) ×160mm	1.4 (17G) ×130mm
SBN-1820C	1.2 (18G) ×200mm	1.4 (17G) ×170mm
SBN-2010C	0.9 (20G) ×100mm	1.1 (19G) ×70mm
SBN-2016C	0.9 (20G) ×160mm	1.1 (19G) ×130mm
SBN-2020C	0.9 (20G) ×200mm	1.1 (19G) ×170mm

※ Apply one reference number (REF) for both the semi-automatic biopsy needle and co-axial biopsy device to simplify ordering and ensure that the co-axial biopsy device is compatible with the semi-automatic biopsy needle.

※ If you have further questions or need other specifications, please contact us.



PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca
www.promisemed-global.com