


JIANGSU HITEC MEDICAL CO., LTD

Hitec Solutions
for Hospital and Home Care
www.hitecare.com




About Us:

Jiangsu Hitec Medical Co., Ltd, located in Nantong city, Jiangsu province is the original manufacturer of Respiratory, Anesthesia, Urological and Infusion products, with the headquarter (Hitec Medical Co., Ltd) located in Shanghai and brand name HITECARE.

Hitecare products have been exporting to 125 customers over 65 countries in Europe, South America, Middles East, Asia, Africa and Australia. Hitecare reputation has been well established in the world market.

Jiangsu Hitec Medical Co.,Ltd, covers an area of 21,307.58 square meters, with 5,000 square meters of Class 100,000 cleanrooms and a total investment of 120 million RMB.



All products are **MDR CE** certified and comply with **ISO 13485** standards registered with TÜV SÜD Product Service GmbH (certification body, identification number 0123).

The Hitecare team is committed to ensuring that every product performs at its best to meet the needs of hospital and home care.

Email: sheng@hitecmed.com

Mobile: 00 86 13472533440

Contact: Jackie Sheng

Anesthesia

01

Standard Endotracheal tube	01
Oral preformed Endotracheal tube	02
Nasal preformed Endotracheal tube	02
Reinforced endotracheal tube	03
Endotracheal tube with suction lumen	03
Intubating stylet	04
Endobronchial tube	04
Endobronchial blocker tube	05
Tracheostomy tube	06
Laryngeal mask	07
Oropharyngeal airway	08
Nasopharyngeal airway	08
Suction catheter	09
Closed suction catheter	09
Suction connecting tube	10

Urology

03

Silicone foley catheter	22
Latex foley catheter	23
Latex foley catheter	24
Latex foley catheter	25
Hydrophilic coated Intermittent catheter	26
Nelaton catheter	27

Respiratory

02

Oxygen mask	11
Nebulizer	12
Non-rebreathing oxygen mask	13
Multi-vent mask /Exchanging type	14
Multi-vent mask /Rotation type	15
Nasal Oxygen cannula	16
Catheter Mount	17
Breathing circuit	17
Air cushion face mask	18
Anesthesia face mask Silicone	18
CPAP Mask	19
Heat and Moisture exchange filter	20
Bacteria filter	20
Manual Resuscitator	21

Infusion

04

Urine bag	28
Leg bag	29
Infusion and Transfusion set	30
Blood transfusion set	31




Standard Endotracheal tube



Product	Size ID(mm)	Ref.code	Product	Size ID(mm)	Ref.code
Uncuffed	2.0	B0101	Cuffed	2.0	B0201
	2.5	B0102		2.5	B0202
	3.0	B0103		3.0	B0203
	3.5	B0104		3.5	B0204
	4.0	B0105		4.0	B0205
	4.5	B0106		4.5	B0206
	5.0	B0107		5.0	B0207
	5.5	B0108		5.5	B0208
	6.0	B0109		6.0	B0209
	6.5	B0110		6.5	B0210
	7.0	B0111		7.0	B0211
	7.5	B0112		7.5	B0212
	8.0	B0113		8.0	B0213
	8.5	B0114		8.5	B0214
	9.0	B0115		9.0	B0215
	9.5	B0116		9.5	B0216

Benefits

- Bevel tip: rounded tip ensures atraumatic and smooth intubation and facilitates fiberoptic intubation.
- Murphy eye: smoothly formed to avoid vocal damage during intubation.
- Soft intubation: A soft transition between cuff and shaft without sharp edges provides atraumatic and smooth in- and extubation.
- Effective seal: cylindrical shaped and thin walled High volume Low pressure(HVLP) cuff for soft and effective sealing.
- Depth marking: double depth marking helps to guide intubation by laryngoscopy ensuring accurate placement.
- Radiopaque x-ray line: facilitates verification of tube positioning visible to the tip.
- Accurate position control: grading in 1cm steps provides reliable control of the tube position.
- 15mm connector, assures safety connection with filters and circuits.
- Medical grade PVC, soft, thermosensitive, kink-resistant ensures optimal anatomical adaptation.





Oral preformed Endotracheal tube



Product	Size ID(mm)	Ref.code	Product	Size ID(mm)	Ref.code
Oral preformed Uncuffed	3.0	B0303	Oral preformed Cuffed	3.0	B0403
	3.5	B0304		3.5	B0404
	4.0	B0305		4.0	B0405
	4.5	B0306		4.5	B0406
	5.0	B0307		5.0	B0407
	5.5	B0308		5.5	B0408
	6.0	B0309		6.0	B0409
	6.5	B0310		6.5	B0410
	7.0	B0311		7.0	B0411
	7.5	B0312		7.5	B0412
	8.0	B0313		8.0	B0413
	8.5	B0314		8.5	B0414
	9.0	B0315		9.0	B0415
	9.5	B0316		9.5	B0416


Nasal preformed Endotracheal tube



Product	Size ID(mm)	Ref.code	Product	Size ID(mm)	Ref.code
Nasal preformed Uncuffed	3.0	B0503	Nasal preformed Cuffed	3.0	B0603
	3.5	B0504		3.5	B0604
	4.0	B0505		4.0	B0605
	4.5	B0506		4.5	B0606
	5.0	B0507		5.0	B0607
	5.5	B0508		5.5	B0608
	6.0	B0509		6.0	B0609
	6.5	B0510		6.5	B0610
	7.0	B0511		7.0	B0611
	7.5	B0512		7.5	B0612
	8.0	B0513		8.0	B0613
	8.5	B0514		8.5	B0614
	9.0	B0515		9.0	B0615
	9.5	B0516		9.5	B0616

Benefits

- The preformed bend lies comfortably at the nares and can be cut shorter for easier suction.
- The design ensures good access to the face allowing the surgeon to observe the effect of the maxilla manipulation on facial symmetry.





Reinforced endotracheal tube



Product	Size ID(mm)	Ref.code	Product	Size ID(mm)	Ref.code
Reinforced Uncuffed	3.0	B0703	Reinforced Cuffed	3.0	B0803
	3.5	B0704		3.5	B0804
	4.0	B0705		4.0	B0805
	4.5	B0706		4.5	B0806
	5.0	B0707		5.0	B0807
	5.5	B0708		5.5	B0808
	6.0	B0709		6.0	B0809
	6.5	B0710		6.5	B0810
	7.0	B0711		7.0	B0811
	7.5	B0712		7.5	B0812
	8.0	B0713		8.0	B0813
	8.5	B0714		8.5	B0814
	9.0	B0715		9.0	B0815
	9.5	B0716		9.5	B0816

Benefits

- High strength compression resistant spring, no matter how the patient's position is changed, the tube will not be deflated or distorted.
- Facilitates some special positions for surgery, prone or back surgery.
- Provides the perfect balance of kink resistance, flexibility and stiffness.


Endotracheal tube with suction lumen

Product	Size ID(mm)	Ref.code
with suction lumen	6.0	B0909
	6.5	B0910
	7.0	B0911
	7.5	B0912
	8.0	B0913
	8.5	B0914
	9.0	B0915
	9.5	B0916



Benefits

- A suction hold above the cuff, allows continuously suctioning of subglottic secretions, which will greatly reduce the risk of ventilator associated pneumonia(VAP), consequently reduce patient risk and medical-care cost.



Intubating stylet



Product	Size(Fr)	Ref.code
Intubating stylet	06	DS002
	10	DS004
	14	DS023

Bougie


Product	Size	Type	Ref.code
Bougie	6	Solid	Z0101
	10	Solid	Z0102
	15	Solid	Z0103
	10	Hollow	Z0202
	15	Hollow	Z0203



Benefits

- The outer casing is made of PVC.
- Easy to insert and remove tubes with smooth surface.
- Rounded head tip to minimize damage to the tracheal wall.

Endobronchial tube



Product	Size ID(Fr)	Ref.code	Product	Size ID(Fr)	Ref.code
Left	28	V01010105	Right	28	V01010205
	32	V01010101		32	V01010201
	35	V01010102		35	V01010202
	37	V01010103		37	V01010203
	39	V01010104		39	V01010204
	41	V01010106		41	V01010206

Benefits

- The bronchial cuff performs better for sealing the tracheal secretions.
- Pilot balloons fitted on the inflation tube to indicate inflation of the bronchial cuff and tracheal cuff.
- Connectors assure safety connection with filters and circuits.
- For left or right lung.





Endobronchial blocker tube



Product	Size ID(FR)	Ref.code
Blocker	05	C-I01
	07	C-I02
	09	C-I03

Benefits

- Endobronchial blocker tube is a device which can be inserted down a tracheal tube after tracheal intubation so as to block off the right or left main bronchus of the lungs in order to be able to achieve a controlled one-sided ventilation of the lungs in thoracic surgery.
- The stability of cuff balloon's internal pressure is better.
- The connector enables the connection of an endobronchial blocker tube to various tracheal tube or breathing circuits.





Tracheostomy tube



Product	Size ID(mm)	Ref.code	Product	Size ID(mm)	Ref.code
Cuffed	3.0	L0103	Uncuffed	3.0	L0203
	3.5	L0104		3.5	L0204
	4.0	L0105		4.0	L0205
	4.5	L0106		4.5	L0206
	5.0	L0107		5.0	L0207
	5.5	L0108		5.5	L0208
	6.0	L0109		6.0	L0209
	6.5	L0110		6.5	L0210
	7.0	L0111		7.0	L0211
	7.5	L0112		7.5	L0212
	8.0	L0113		8.0	L0213
	8.5	L0114		8.5	L0214
	9.0	L0115		9.0	L0215
	9.5	L0116		9.5	L0216

Benefits

- Belt: patients with abrasions or burns on the neck are less likely to adhere to the trauma.
- Pilot balloon: balloon fitted to an inflating tube to indicate inflation of the cuff.
- Main tube: designed with light weight medical grade PVC.
- Cuff: cylindrical cuff allows for a large contact area, low pressure, minimal damage and good sealing.
- Full length Radio-opaque line for correct placement detection.
- Smooth rounded tip of the Obturator reduces trauma during insertion.
- Printed neck plate with size information for easy identification.



Laryngeal mask



Product	Size ID	Ref.code	Product	Size ID	Ref.code
Single use PVC	1.0 #	O0101	Single use PVC Reinforced	1.0 #	O0401
	1.5 #	O0102		1.5 #	O0402
	2.0 #	O0103		2.0 #	O0403
	2.5 #	O0104		2.5 #	O0404
	3.0 #	O0105		3.0 #	O0405
	4.0 #	O0106		4.0 #	O0406
	5.0 #	O0107		5.0 #	O0407
Single use Silicone	1.0 #	O0201	Single use Silicone Reinforced	1.0 #	O0301
	1.5 #	O0202		1.5 #	O0302
	2.0 #	O0203		2.0 #	O0303
	2.5 #	O0204		2.5 #	O0304
	3.0 #	O0205		3.0 #	O0305
	4.0 #	O0206		4.0 #	O0306
	5.0 #	O0207		5.0 #	O0307
Reusable Silicone	1.0 #	O0201C	Reusable Silicone Reinforced	1.0 #	O0301C
	1.5 #	O0202C		1.5 #	O0302C
	2.0 #	O0203C		2.0 #	O0303C
	2.5 #	O0204C		2.5 #	O0304C
	3.0 #	O0205C		3.0 #	O0305C
	4.0 #	O0206C		4.0 #	O0306C
	5.0 #	O0207C		5.0 #	O0307C

Benefits

- Mask non-displacement, easy insertion and prevention of structural damage to human tissue.
- Airway tube designed in lightweight medical grade PVC or Silicone for patient comfort.
- One-way valve for inflation and deflation.
- Easy to use, rapid and easier to maintain the airway.
- No laryngoscopy is required, making LMA placement less difficult and more successful for the beginner.
- Avoiding endotracheal mucosal damage.
- The airway is established for autonomous and controlled ventilation.





Oropharyngeal airway



Product	Size ID	Ref-code	Color coding
Guedel type	40 mm 000 #	Q0101	Pink
	50 mm 00 #	Q0102	Dark blue
	60 mm 0 #	Q0103	Black
	70 mm 1 #	Q0104	White
	80 mm 2 #	Q0105	Green
	90 mm 3 #	Q0106	Yellow
	100 mm 4 #	Q0107	Red
	110 mm 5 #	Q0108	Light blue
	120 mm 6 #	Q0109	Orange

Benefits

- Atraumatic softly rounded tip ensures atraumatic and smooth intubation.
- Ergonomic design of the curved tube assists intubation and enhances patient comfort.
- The flange keeps the oropharyngeal airway in place and does not slip easily.
- Color coded for easy size identification, prevents oropharyngeal airway from being damaged by teeth.
- Prevents the tongue from blocking the upper airway.


Nasopharyngeal airway



Product	Size ID(mm)	Ref.code	Product	Size ID(mm)	Ref.code	Product	Size ID(mm)	Ref.code
Silicone Green	2.5	R0401LG	PVC Transparent Green	2.5	R0101	PVC Blue	2.5	R0101B
	3.0	R0402LG		3.0	R0102		3.0	R0102B
	3.5	R0403LG		3.5	R0103		3.5	R0103B
	4.0	R0404LG		4.0	R0104		4.0	R0104B
	4.5	R0405LG		4.5	R0105		4.5	R0105B
	5.0	R0406LG		5.0	R0106		5.0	R0106B
	5.5	R0407LG		5.5	R0107		5.5	R0107B
	6.0	R0408LG		6.0	R0108		6.0	R0108B
	6.5	R0409LG		6.5	R0109		6.5	R0109B
	7.0	R0410LG		7.0	R0110		7.0	R0110B
	7.5	R0411LG		7.5	R0111		7.5	R0111B
	8.0	R0412LG		8.0	R0112		8.0	R0112B
	8.5	R0413LG		8.5	R0113		8.5	R0113B
	9.0	R0414LG		9.0	R0114		9.0	R0114B

Benefits

- Smooth bevel tip, less trauma, reduce the friction on the nasal tissue when intubation tip and passes through the nasal cavity smoothly.
- Flared end, softly round, helps keep the nasopharyngeal airway in place and does not slip easily.
- With size identification on tubing, easy for choosing.
- The nasopharyngeal airway does not cause nausea and vomiting, and does not affect speaking, opening, closing and swallowing.
- Suitable for minor oral surgery.
- Stomach tube can be placed through the nasopharyngeal airway and is generally not easy to curl.





Suction catheter

Plain type




T type



Cap-cone type



Y type



Size(FR/CH)	Color coding		Plain type	T type	Cap-cone	Y type
06	Light green	Ref.code	K0101	K0201	K0301	K0401
08	Blue		K0102	K0202	K0302	K0402
10	Black		K0103	K0203	K0303	K0403
12	White		K0104	K0204	K0304	K0404
14	Green		K0105	K0205	K0305	K0405
16	Orange		K0106	K0206	K0306	K0406
18	Red		K0107	K0207	K0307	K0407
20	Yellow		K0108	K0208	K0308	K0408

Benefits

- Catheter with X-ray line(optional), easy for checking position of catheter, flexible maintains efficient suction with ease of insertion and withdrawal.
- The connector is effective at presenting contact with fluids, thus helping to minimizing contamination.
- Staggered and opposite eyes reduce potential for invagination or blocked suction.
- Smoothly rounded tip for less trauma.
- Color coded connector for size identification.

Closed suction catheter

Benefits

- Connector: Double swivel connectors reduce strain on ventilator tubing.
- Medicine injection port: The airway can be continuously moistened by drug injection to reduce the viscosity of sputum and reduce the incidence of sputum blockage in the artificial airway.
- Irrigation port: It is well sealed and with anti-reflux design, which avoids extravasation when flushing .
- Vacuum control: Rotating manifold locks into place with a click for reassurance that circuit remain sealed.
- Isolation valve: Ergonomic suction control switch with security lock and an easy single hand operation design.
- Sleeve: Tactile sleeve for catheter advancement and prevents cross infection.
- Suction catheter: The catheter with X-ray line, easy for checking position of catheter and depth markings for easy and proper depth suctioning.




Duration	Size(FR)	Ref.code	Color coding
24 Hours	06	X01010101	Red
	08	X01010201	Light blue
	10	X01010301	Black
	12	X01010401	White
	14	X01010501	Green
	16	X01010601	Orange
	18	X01010701	Red
	20	X01010801	Yellow
72 Hours	06	X01010103	Red
	08	X01010203	Light blue
	10	X01010303	Black
	12	X01010403	White
	14	X01010503	Green
	16	X01010603	Orange
	18	X01010703	Red
	20	X01010803	Yellow



Suction connecting tube

Product	Type	Vacuum control	Ref.code
Yankauer Handle	Plain tip	without	C020002
	Plain tip	with	C020005
	Crown tip	without	C020003
	Crown tip	with	C020006



Product	ID(mm)	OD(mm)	CH		1.8m	2.0m	2.5m	3m	3.6m
Suction connecting tube	5.5	8.2	22	Ref. code	C010100	C010700	C011300	C011900	C012500
	6.0	9.0	24		C010200	C010800	C011400	C012000	C012600
	7.0	10.0	28		C010400	C011000	C011600	C012200	C012800



Standard female connector



Medium flared connector



One cut to fit connector



Transition connector without cap




Transition connector with cap

Length	ID(mm)	OD(mm)	CH		Plain tip without vacuum control	Crown tip without vacuum control	Plain tip with vacuum control	Crown tip with vacuum control
1.8m	5.5	8.2	22	Ref. code	C030102	C030103	C030105	C030106
	6.0	9.0	24		C030202	C030203	C030205	C030206
	7.0	10.0	28		C030402	C030403	C030405	C030406
2.0m	5.5	8.2	22		C030702	C030703	C030705	C030706
	6.0	9.0	24		C030802	C030803	C030805	C030806
	7.0	10.0	28		C031002	C031003	C031005	C031006
2.5m	5.5	8.2	22		C031302	C031303	C031305	C031306
	6.0	9.0	24		C031402	C031403	C031405	C031406
	7.0	10.0	28		C031602	C031603	C031605	C031606
3.0m	5.5	8.2	22		C031902	C031903	C031905	C031906
	6.0	9.0	24		C032002	C032003	C032005	C032006
	7.0	10.0	28		C032202	C032203	C032205	C032206

Benefits

- Handle: The forming angle of handle is designed to make it more convenient and humanized in clinic use and hex-arris design to avoid blocking the tube due to high pressure.
- Connecting tube with optimal draping qualities for maintaining a wall thickness to prevents the tube from collapsing when the tube is used under high negative pressure.
- Tip with a large opening surrounded by a bulbous head and is designed to allow effective suction without damaging surrounding tissue.





Benefits

- Mask with ergonomic design, different sizes of mask fit the chins of different ages, seal well.
The mask is highly transparent, which can detect the patient's vomiting or abnormality.
- Adjustable nose clip, easy to fit various nose shapes, and prevent the air flow from damaging the eyes during using.
- Adjustable head belt, fit various head shapes, fix the mask in required position.
- Oxygen tube with star lumen designing, to minimize risk of airflow termination when the tube is kinked.

Product	Color	Length		Pediatric standard(S)	Pediatric elongated(M)	Adult standard(L)	Adult elongated(XL)
Sterile	Green transparent	1.5m	Ref.code	D010101L	D010201L	D010301L	D010401L
		2.0m		D010102L	D010202L	D010302L	D010402L
	White transparent	1.5m		D010101B	D010201B	D010301B	D010401B
		2.0m		D010102B	D010202B	D010302B	D010402B
Non-Sterile	Green transparent	1.5m		D020101L	D020201L	D020301L	D020401L
		2.0m		D020102L	D020202L	D020302L	D020402L
	White transparent	1.5m		D020101B	D020201B	D020301B	D020401B
		2.0m		D020102B	D020202B	D020302B	D020402B



Nebulizer with mask



Benefits

- Mask with ergonomic design, different sizes of mask fit the chins of different ages, seal well. The mask is highly transparent, which can detect the patient's vomiting or abnormality.
- Adjustable nose clip, easy to fit various nose shapes, and prevent the air flow from damaging the eyes during using.
- Adjustable head belt, fit various head shapes, fix the mask in required position.
- Oxygen tube with star lumen designing, to minimize risk of airflow termination when the tube is kinked.
- Nebulizer delivers fast and effective treatment, allowing smaller particles of the medication to be breathed in efficiently.

Product	Color	Length		Pediatric standard(S)	Pediatric elongated(M)	Adult standard(L)	Adult elongated(XL)
Sterile	Green transparent	1.5m	Ref.code	F010101L	F010201L	F010301L	F010401L
		2.0m		F010102L	F010202L	F010302L	F010402L
	White transparent	1.5m		F010101B	F010201B	F010301B	F010401B
		2.0m		F010102B	F010202B	F010302B	F010402B
Non-Sterile	Green transparent	1.5m		F020101L	F020201L	F020301L	F020401L
		2.0m		F020102L	F020202L	F020302L	F020402L
	White transparent	1.5m		F020101B	F020201B	F020301B	F020401B
		2.0m		F020102B	F020202B	F020302B	F020402B

Nebulizer with mouth piece

Product	Color	Length	Ref.code
Sterile	Green transparent	1.5m	F030301L
		2.0m	F030302L
	White transparent	1.5m	F030301B
		2.0m	F030302B
Non-Sterile	Green transparent	1.5m	F040301L
		2.0m	F040302L
	White transparent	1.5m	F040301B
		2.0m	F040302B



Non-rebreathing oxygen mask



Benefits

- Mask with ergonomic design, different sizes of mask fit the chins of different ages, seal well.
The mask is highly transparent, which can detect the patient's vomiting or abnormality.
- Adjustable nose clip, easy to fit various nose shapes, and prevent the air flow from damaging the eyes during using.
- Adjustable head belt, fit various head shapes, fix the mask in required position.
- Oxygen tube with star lumen designing, to minimize risk of airflow termination when the tube is kinked.
- Three silicone valves are used to prevent rebreathing of exhaled CO₂.
- Reservoir bag provide high oxygen concentration, being up to 90% and even higher, for better therapy effect.
Reservoir bag available with 1000ml and 750ml.

Product	Color	Length		Pediatric standard(S)	Pediatric elongated(M)	Adult standard(L)	Adult elongated(XL)
Sterile	Green transparent	1.5m	Ref.code	E010101L	E010201L	E010301L	E010401L
		2.0m		E010102L	E010202L	E010302L	E010402L
	White transparent	1.5m		E010101B	E010201B	E010301B	E010401B
		2.0m		E010102B	E010202B	E010302B	E010402B
Non-Sterile	Green transparent	1.5m		E020101L	E020201L	E020301L	E020401L
		2.0m		E020102L	E020202L	E020302L	E020402L
	White transparent	1.5m		E020101B	E020201B	E020301B	E020401B
		2.0m		E020102B	E020202B	E020302B	E020402B



Multi-vent mask /Exchanging type



Benefits

- Mask with ergonomic design, different sizes of mask fit the chins of different ages, seal well. The mask is highly transparent, which can detect the patient's vomiting or abnormality.
- Adjustable nose clip, easy to fit various nose shapes, and prevent the air flow from damaging the eyes during using.
- Adjustable head belt, fit various head shapes, fix the mask in required position.
- Oxygen tube with star lumen designing, to minimize risk of airflow termination when the tube is kinked.
- The humidifier cup is connected to the humidifier at one end to make the dry oxygen produced by the oxygen concentrator more moist and comfortable for the patient to breathe, and its hollow interior allows for the replacement of multiple oxygen concentration control connectors to provide the appropriate oxygen concentration for patients with different needs (Blue 24%, Yellow 28%, White 31%, Green 35%, Pink 40%, Orange 50%, Red 60%).

Product	Length	Color		Pediatric standard(S)	Pediatric elongated(M)	Adult standard(L)	Adult elongated(XL)
Sterile	2.0m	Green transparent	Ref.code	G010102L	G010202L	G010302L	G010402L
		White transparent		G010102B	G010202B	G010302B	G010402B
Non-Sterile	2.0m	Green transparent		G020102L	G020202L	G020302L	G020402L
		White transparent		G020102B	G020202B	G020302B	G020402B



Multi-vent mask /Rotation type




Benefits

- Mask with ergonomic design, different sizes of mask fit the chins of different ages, seal well.
The mask is highly transparent, which can detect the patient's vomiting or abnormality.
- Adjustable nose clip, easy to fit various nose shapes, and prevent the air flow from damaging the eyes during using.
- Adjustable head belt, fit various head shapes, fix the mask in required position.
- Oxygen tube with star lumen designing, to minimize risk of airflow termination when the tube is kinked.
- Multi-vent set: connecting a humidifier makes the patient's breathing more uncomfortable because the output of the oxygen generator is quite dry (Green connector for concentration 24%, 26%, 28% and 30%, White for 35%, 40%, 50%).

Product	Length	Color		Pediatric standard(S)	Pediatric elongated(M)	Adult standard(L)	Adult elongated(XL)
Sterile	2.0m	Green transparent	Ref.code	G030102L	G030202L	G030302L	G030402L
		White transparent		G030102B	G030202B	G030302B	G030402B
Non-Sterile	2.0m	Green transparent		G040102L	G040202L	G040302L	G040402L
		White transparent		G040102B	G040202B	G040302B	G040402B



Nasal Oxygen cannula



1. Air supply connector
2. Tube B
3. Tee joint
4. Adjusting ring
5. Tube A
6. Nostril(PVC,silicone)

Benefits


- Connector of oxygen tubing can be applied to various oxygen supply equipment ends.
- Double channels are more conducive to deliver supplemental oxygen to a patient or person in need of extra oxygen.
- Star lumen designing of oxygen tubing, to minimize risk of airflow termination when the tube is kinked.
- Nostril sucker, soft and smooth surface without sharp edge and object, which does not cause irritation to the patient's nostrils.

PVC nostril sucker (Injection molding)



Product	Color	Length		Tube A 2.7mm / Tube B 5.0mm		Tube A 3.3mm / Tube B 5.0mm	
				Neonate	Infant	Child	Adult
Sterile	Green transparent	1.0m	Ref.code	A010101L	A010301L	A010501L	A010701L
		1.5m		A010102L	A010302L	A010502L	A010702L
		2.0m		A010103L	A010303L	A010503L	A010703L
	White transparent	1.0m		A010101B	A010301B	A010501B	A010701B
		1.5m		A010102B	A010302B	A010502B	A010702B
		2.0m		A010103B	A010303B	A010503B	A010703B
Non-Sterile	Green transparent	1.0m		A020101L	A020301L	A020501L	A020701L
		1.5m		A020102L	A020302L	A020502L	A020702L
		2.0m		A020103L	A020303L	A020503L	A020703L
	White transparent	1.0m		A020101B	A020301B	A020501B	A020701B
		1.5m		A020102B	A020302B	A020502B	A020702B
		2.0m		A020103B	A020303B	A020503B	A020703B

PVC nostril sucker (Dipping molding)



Product	Color	Length		Tube A 2.7mm / Tube B 5.0mm		Tube A 3.3mm / Tube B 5.0mm	
				Neonate	Infant	Child	Adult
Sterile	White transparent	1.0m	Ref.code	A010201B	A010401B	A010601B	A010801B
		1.5m		A010202B	A010402B	A010602B	A010802B
		2.0m		A010203B	A010403B	A010603B	A010803B
Non-Sterile	White transparent	1.0m		A020201B	A020401B	A020601B	A020801B
		1.5m		A020202B	A020402B	A020602B	A020802B
		2.0m		A020203B	A020403B	A020603B	A020803B

Catheter Mount



Product	Type	Ref.code
Catheter Mount	collapsible	V1201
	smoothbore	V1001
	corrugated	V1101




Breathing circuit

Product	Type	Ref.code	Product	Type	Ref.code
Child	1.6m corrugated tube with 1L gas bag	V04030200	Adult	1.6m corrugated tube with 2L gas bag	V04040100
	1.6m Flexible tube with 1L gas bag	V06030200		1.6m Flexible tube with 2L gas bag	V06040100
	1.6m smooth bore tube with 1L gas bag	V05030200		1.6m smooth bore tube with 2L gas bag	V05040100
	1.6m corrugated tube with 2 water traps, with 0.5m extra limb	V07030000		1.6m corrugated tube with 2 water traps, with 0.5m extra limb	V07040000
	1.6m Flexible tube with 2 water traps, with 0.5m extra limb	V09030000		1.6m Flexible tube with 2 water traps, with 0.5m extra limb	V09040000
	1.6m smooth bore tube with 2 water traps, with 0.5m extra limb	V08030000		1.6m smooth bore tube with 2 water traps, with 0.5m extra limb	V08040000

Benefits

- Conservation of Gas and Heat: Closed systems particularly conserve anesthetic gases and patient body heat, making procedures more efficient and comfortable for the patient.
- Safety: Modern circuits come equipped with fail-safes and monitors that reduce the risk of delivering a hypoxic mix (low oxygen content) to the patient.
- Flexibility: The design of these circuits allows for easy modification to cater to individual patient needs, be it adults, children, or infants.
- Humidification: The circuit can provide moisture to the gases the patient breathes, reducing the risk of a dry respiratory tract.
- CO2 Removal: Efficiently removes carbon dioxide from the patient's exhaled gases, ensuring a fresh breath in the next cycle.





Air cushion face mask

Benefits

- Fixed hook ring, color coded for easy size identification and fixation.
- The mask cup is highly transparent, which for the observation of natural breathing blood stained and patient's vomit. It also can be connected with multiply medical equipment.
- Inflation valve with good air tightness to make it possible to inflate the air with cushion.
- Use air cushion type design of joint face, soft and comfortable, high air tightness, avoid air leakage at the start of the clinical procedure.



Product	Type		Neonate	Infant	Child	Small adult	Medium adult	Large adult	Ultra Large adult
Air cushion face mask	Horizontal	Ref.code	T0201	T0202	T0203	T0204	T0205	T0206	T0207
	Vertical		T0101	T0102	T0103	T0104	T0105	T0106	T0107

Anesthesia face mask Silicone


Benefits

- Latex Free, PVC Free.
- Autoclavable for reuse.



Product	Size	Ref.code
Air cushion face mask	0# Neonate	T0301
	1# Infant	T0302
	2# Child	T0303
	3# Small adult	T0304
	4# Medium adult	T0305
	5# Large adult	T0306
	6# Ultra Large adult	T0307





CPAP Mask




Benefits

- Deliver continuous positive airway pressure to patient, thereby enabling the patient to inhale fresh air with each respiration.
- The hollow frame design is easier to clean and maintain.
- Anti-sliding headgear with ideal elasticity and air permeability to fulfill fantastic comfort, suitable for sensitive skins.
- The mask is highly transparent, which for the observation of natural breathing bloodstained and patients' vomit. It is also soft and comfortable, high air tightness, avoid air leakage.
- It can be connected with multiple medical equipment to provide patients with oxygen.
- DEHP FREE, LATEX FREE.

Product	CPAP MASK	
Size	M	L
Ref.code	S020102	S030102
Silicone mask dimension ID	77.3*55 mm	88.3*63.5 mm
Mask cover dimension	200*119*107 mm	
Weight	133g	140g
Mask connector	22mm conical	
Positive pressure	4~30cm H ₂ O	
Air flow resistance	Pressure against valve opening≤ 2cm H ₂ O	
Anti-suffocation valve open	Pressure against valve closing≤ 2.5cm H ₂ O	
Anti-suffocation valve opening and closing pressure	≤ 43dB	
Working environment temperature	5℃ ~ 40℃	
Working environment humidity	Relative humidity 10%~95%, no condensation	
Storage temperature	-20℃ ~ 40℃	
Storage humidity	Relative humidity 10%~95%, no condensation	





Heat and Moisture exchange filter



Benefits

- Providing optimal moisture and temperature output with low resistance to flow and bi-directional filtration with bacterial/viral efficiency offering cross-contamination protection for patients when the clinical gas pass through.
- The standard 15M/22F machine end allows for better connection to circuits
- Patient end 15F/22M.
- With Gas sampling port to meet clinical medical needs.
- Filter core consists of a filter membrane and a wetting paper made of a double layer of corrugated paper. This filter membrane effectively reduces the spread of microorganisms and other particles in the respiratory system. The paper adds sufficient wet and heat to the patient's inhaled air without causing resistance to gas delivery, improving the patient's lung function.

Product		Type	Infant	Child	Adult
Heat and Moisture exchange filter	Sterile	ABS	RB059	RB023	RB020
		PP	RB060	RB024	RB021


Bacteria filter



Product		Type	Adult
Bacteria filter	Sterile	ABS	GL020
		PP	GL021

Benefits

- Used for bi-directional filtration with bacterial/viral efficiency offering cross-contamination protection for patients when the clinical gas pass through.
- The standard 15M/22F machine end allows for better connection to circuits.
- Patient end 15F/22M.
- With Gas sampling port to meet clinical medical needs.
- Filter core consists of a filter membrane which effectively reduces the spread of microorganisms and other particles in the respiratory system, without causing resistance to gas delivery, improving the patient's lung function.





Manual Resuscitator



Benefits

- Simple structure, easy to carry and good ventilation.
- Safety valve to prevent rupture of alveoli due to excessive pressure, ensuring patient safety.
- PVC manual resuscitator is for single use.
- Silicone manual resuscitator main parts are autoclavable can be reused for 20 times.

Clean and Re-sterilize:

- Wash the Silicone resuscitator with 0.2% Javel water and rinse in clean water.
- Autoclave the Silicone mask, Silicone bag, valve for 5 minutes at 134 °C or 15 minutes at 121 °C.
- The reservoir bag, oxygen tubing, mouth opener(if provided), Guedel airway(if provided) are single use, can not be washed or autoclaved.

Product		Manual Resuscitator		
Size		Infant	Child	Adult
Ref.code	PVC	Z0867	Z0717	Z1051
	Silicone	Z0531	Z0109	Z0132
Body weight		≤ 10 kg	10kg~40kg	> 40 kg
Stroke volume		150ml	400ml	800ml
Resuscitator volume		280ml	600ml	1650ml
Dimensions(Length x Diameter)		225×85mm	262×105mm	325×132mm
Pressure limiting valve		40 cmH ₂ O	40 cmH ₂ O	60 cmH ₂ O
Dead space(patient valve)		7 ml	7 ml	7 ml
Inspiratory resistance		<5cm H ₂ O(at 50L/min)		
Expiratory resistance		<5cm H ₂ O(at 50L/min)		
Bag reservoir volume		1600 ml	1600 ml	2000 ml
Patient connector		ISO5356-1: Ø 22/15mm		
Recommended operating condition		-18 °C ~ +50 °C		
Storage		-40 °C ~ +60 °C		
Supplementary oxygen and delivered oxygen concentration(≥85%)		V _T :60ml Frequency:20BPM Flow rate:15L/min	V _T :150ml Frequency:25BPM Flow rate:15L/min	V _T :600ml Frequency:12BPM Flow rate:15L/min


Silicone foley catheter



Product	Size(FR)	Balloon capacity	Ref.code	Color coding	Product	Size(FR)	Balloon capacity	Ref.code	Color coding
2 way Straight tip	06	3 ml	I010601	Green	2 way Tieman tip	12	5-10 ml	I021202	White
	08	3-5 ml	I010801	Blue		14	5-10 ml	I021402	Green
	10	3-5 ml	I011001	Black		16	5-10 ml	I021602	Orange
	12	5-10 ml	I011202	White		16	30 ml	I021605	Orange
	14	5-10 ml	I011402	Green		18	30 ml	I021805	Red
	16	5-10 ml	I011602	Orange		20	30 ml	I022005	Yellow
	16	30 ml	I011605	Orange		22	30 ml	I022205	Purple
	18	30 ml	I011805	Red		24	30 ml	I022405	Blue
	20	30 ml	I012005	Yellow		26	30 ml	I022605	Pink
	22	30 ml	I012205	Purple					
	24	30 ml	I012405	Blue					
	26	30 ml	I012605	Pink					
3 way Straight tip	14	5-10 ml	I031402	Green	3 way Tieman tip	14	5-10 ml	I041402	Green
	16	5-10 ml	I031602	Orange		16	5-10 ml	I041602	Orange
	16	30 ml	I031605	Orange		16	30 ml	I041605	Orange
	18	30 ml	I011805	Red		18	30 ml	I041805	Red
	20	30 ml	I031805	Yellow		20	30 ml	I042005	Yellow
	22	30 ml	I032205	Purple		22	30 ml	I042205	Purple
	24	30 ml	I032405	Blue		24	30 ml	I042405	Blue
	26	30 ml	I032605	Pink		26	30 ml	I042605	Pink

Benefits

- The shaft's inner layer strength aids in insertion by providing smooth, uniform walls and resists collapse during aspiration. The shaft with X-ray line, easy for checking position of catheter
- Drainage funnel round with optimal internal diameter which resists obstruction and potential for encrustation and provides secure connection to most drainage tubes.
- Inflation access port is designed ergonomically which enhance speed of inflation and deflation.
- The balloon with correct anatomical shape to hold tip in proper upright position more permanently, which is designed to retain its shape with good concentricity.
- Tip with X-ray, easy for checking position of catheter tip, which is smooth rounded so as to reduce trauma during insertion.
- The valve accommodates most luer slip and luer lock syringes.





Latex foley catheter

Features

Natural latex, medical-grade

Silicone coated

Available with balloon in variety of capacity



2 Way Pediatric Rubber valve,length 270mm

Size(FR/CH)	Color coding	Ref.code	3ml	3-5ml
06	Pink		H030601	
08	Black		H030801	H030802
10	Gray		H031001	H031002



2 Way Standard Rubber valve,length 400mm

Size(FR/CH)	Color coding	Ref.code	5-10ml	5-15ml	30ml
11	white		H031204	H031205	H031207
11	Green		H031404	H031405	H031407
11	Orange		H031604	H031605	H031607
11	Red		H031804	H031805	H031807
21	Yellow		H032004	H032005	H032007
21	Purple		H032204	H032205	H032207
21	blue			H032405	H032407
21	pink			H032605	H032607
21	Brown			H032805	H032807
31	Gray			H033005	H033007

3 Way Standard Rubber valve,length 400mm



Size(FR/CH)	Ref.code	30ml
16		H091607
18		H091807
20		H092007
22		H092207
24		H092407
26		H092607
28		H092807





Latex foley catheter



2 Way Pediatric Plastic valve,length 270mm

Size(FR/CH)	Color coding	Ref.code	3ml	3-5ml
06	Pink		H040601	
08	Black		H040801	H040802
10	Gray		H041001	H041002



2 Way Standard, Plastic valve,length 400mm

Size(FR/CH)	Color coding	Ref.code	5-10ml	5-15ml	30ml
12	white		H041204	H041205	H041207
14	Green		H041404	H041405	H041407
16	Orange		H041604	H041605	H041607
18	Red		H041804	H041805	H041807
20	Yellow		H042004	H042005	H042007
22	Purple		H042204	H042205	H042207
24	blue			H042405	H042407
26	pink			H042605	H042607
28	Brown			H042805	H042807
30	Gray			H043005	H043007



3 Way Standard, Plastic valve,length 400mm

Size(FR/CH)	Color coding	Ref.code	30ml
16	Orange		H101607
18	Red		H101807
20	Yellow		H102007
22	Purple		H102207
24	blue		H102407
26	pink		H102607
28	Brown		H102807

2 Way Female, Rubber valve,length 240mm




Size(FR/CH)	Color coding	Ref.code	5-10ml	30ml
12	white		H051204	H051207
14	Green		H051404	H051407
16	Orange		H051604	H051607
18	Red		H051804	H051807
20	Yellow		H052004	H052007
22	Purple		H052204	H052207
24	blue			H052407

2 Way Female, Plastic valve,length 240mm



Size(FR/CH)	Color coding	Ref.code	5-10ml	30ml
12	white		H061204	H061207
14	Green		H061404	H061407
16	Orange		H061604	H061607
18	Red		H061804	H061807
20	Yellow		H062004	H062007
22	Purple		H062204	H062207
24	blue			H062407





Latex foley catheter

2 Way Tiemann, Rubber valve,length 400mm

Size(FR/CH)	Color coding		5-10ml	30ml
12	white	Ref.code	H071204	H071207
14	Green		H071404	H071407
16	Orange		H071604	H071607
18	Red		H071804	H071807
20	Yellow		H072004	H072007
22	Purple		H072204	H072207
24	blue			H072407



2 Way Tiemann, Plastic valve,length 400mm

Size(FR/CH)	Color coding		5-10ml	30ml
12	white	Ref.code	H081204	H081207
14	Green		H081404	H081407
16	Orange		H081604	H081607
18	Red		H081804	H081807
20	Yellow		H082004	H082007
22	Purple		H082204	H082207
24	blue			H082407

3 Way Tiemann, Rubber valve,length 400mm

Size(FR/CH)	Color coding		30ml
16	Orange	Ref.code	H121607
18	Red		H121807
20	Yellow		H122007
22	Purple		H122207
24	blue		H122407

3 Way Tiemann, Plastic valve,length 400mm

Size(FR/CH)	Color coding		30ml
16	Orange	Ref.code	H131607
18	Red		H131807
20	Yellow		H132007
22	Purple		H132207
24	blue		H132407

Latex urethral catheter



Size(FR/CH)	Color coding	Ref.code
6	Red	H020600
8		H020800
10		H021000
12		H021200
14		H021400
16		H021600





Hydrophilic coated Intermittent catheter


Gender	Product	Size(FR)	Color coding		Without water sachet	With water sachet	Gender	Type	Size(FR)	Ref.code	Color coding
Male	Straight tip	06	Light green	Ref. code	U010100	U010101	Female (Straight tip)	Without water sachet	06	U011400	Light green
		08	Blue		U010200	U010201			08	U011500	Blue
		10	Black		U010300	U010301			10	U011600	Black
		12	White		U010400	U010401			12	U011700	White
		14	Green		U010500	U010501			14	U011800	Green
		16	Orange		U010600	U010601			16	U011900	Orange
		18	Red		U010700	U010701			18	U012000	Red
		20	Yellow		U010800	U010801			20	U012100	Yellow
	Tiemann tip	06	Light green		U020100	U020101		With water sachet	06	U011401	Light green
		08	Blue		U020200	U020201			08	U011501	Blue
		10	Black		U020300	U020301			10	U011601	Black
		12	White		U020400	U020401			12	U011701	White
		14	Green		U020500	U020501			14	U011801	Green
		16	Orange		U020600	U020601			16	U011901	Orange
		18	Red		U020700	U020701			18	U012001	Red
		20	Yellow		U020800	U020801			20	U012101	Yellow
						Children (Straight tip) (length 250mm)	Without water sachet	06	U012700	Pink	
								08	U012800	Black	
								10	U012900	Gray	
							With water sachet	06	U012701	Pink	
								08	U012801	Black	
								10	U012901	Gray	

Benefits

- Pre-lubricated with the instantly ready-to-use hydrophilic coating, ensure a convenient and simple catheterization every day.
- Designed to reduce friction coefficient and mucosal injury to address comfort due to the super smooth surface.
- Color-coded funnel to easily identify size.
- Not made with natural rubber latex, and phthalate-free.
- The hot-melt punching process, ensures smooth and burr-free drainage apertures.



Nelaton catheter



MALE NELATON CATHETER

NELATON CATHETER TIEMANN TIP

PEDIATRIC NELATON CATHETER



FEMALE NELATON CATHETER

Benefits

- The connector is effective at preventing contact with fluids, thus helping to minimizing contamination and color coded for easy size identification.
- Tip is smooth and rounded so as to reduce trauma during insertion.
- The shaft's inner layer strength aids in insertion by providing smooth, uniform walls and resists collapse.
- Staggered and opposite eyelets reduce potential for invagination or blocked suction.

Gender	Product	Length	Size(FR)	Ref.code	Color coding	Gender	Type	Length	Size(FR)	Ref.code	Color coding
Male	Straight tip	400mm	06	J0101T	Light green	Female	Straight tip	220mm	06	J0114T	Light green
			08	J0102T	Blue				08	J0115T	Blue
			10	J0103T	Black				10	J0116T	Black
			12	J0104T	White				12	J0117T	White
			14	J0105T	Green				14	J0118T	Green
			16	J0106T	Orange				16	J0119T	Orange
			18	J0107T	Red				18	J0120T	Red
			20	J0108T	Yellow				20	J0121T	Yellow
	Tiemann tip	400mm	06	J0201T	Light green				06	J0114T	Light green
			08	J0202T	Blue				08	J0115T	Blue
			10	J0203T	Black				10	J0116T	Black
			12	J0204T	White				12	J0117T	White
			14	J0205T	Green				14	J0118T	Green
			16	J0206T	Orange				16	J0119T	Orange
			18	J0207T	Red				18	J0120T	Red
			20	J0208T	Yellow				20	J0121T	Yellow

Urine bag

Benefits

- With needle free sampling port
- With bed sheet clamp
- With hangers
- T valve with dog house



Economic urine bag



Pediatric urine bag

Product	Size	Ref.code
Economic urine bag	2000ml 90cm Push-pull valve	P01010200
	2000ml 90cm T valve	P01010300
	2000ml 90cm Screw valve	P01010100
Pediatric	100ml	P0701
	200ml	P0702
Luxury Urine bag	2000ml 100cm T valve	P05010300
	4000ml 100cm T valve	P05060300
	5000ml 100cm T valve	P05070300
Simple Luxury Urine Bag	2000ml 90cm T(cross) valve	P02010300
Urine bag with Meter	Bag 2000ml Meter 200ml	P06010301
	Bag 2600ml Meter 200ml	P06050301
	Bag 2600ml Meter 400ml	P06050302
	Bag 2600ml Meter 500ml	P06050303
	Bag 2000ml Meter Two chambers(50+400ml)	P0801



Luxury urine bag



Simple Luxury Urine Bag



Urine bag with Meter



Leg bag

Product	Size	Ref.code
Standard	500ml T valve	P04020300
	750ml T valve	P04030300
	1000ml T valve	P04040300
Non-woven back side	500ml T valve	P04020300W
	750ml T valve	P04030300W
Flocking cloth back side	500ml T valve	P04020300R
	750ml T valve	P04030300R
	1000ml T valve	P04040300R



Standard



**With Non-woven
back side**




**With flocking cloth
back side**



Infusion and Transfusion set

Benefits


- The product can only be used for gravity infusion.
 - If the product is invalid or the package is damaged, it is strictly forbidden to use it.
 - Disposable for use, use immediately after opening, destroy after use.
 - Suitable temperature is 10~40℃.
 - Before use, check whether the package is damaged and whether the sheath has fallen off, otherwise it is not allowed to use.
 - Before use, insert the infusion needle connector tightly to prevent leakage.
- The infusion operation should be implemented and monitored by professional nursing staff.



1. protective cap of closure-piercing device
2. closure-piercing device
3. integral air-inlet with air filter and closure
4. fluid channel
5. drip tube
6. drip chamber
7. fluid filter
8. tubing
9. flow regulator
10. Injection site(latex tube)
11. Needle
12. protective cap of needle
13. Injection site (Y site)

Product	Components	Ref.code	Product	Components	Ref.code
	Spike with air vent	M010601D		Spike with air vent	M020601D
	Drip chamber with filter			Drip chamber with filter	
	Roller			Roller	
	Latex tube injection site			Y injection site	
	With needle			With needle	
	Spike with air vent	M050001D		Spike with air vent	M060001D
	Drip chamber with filter			Drip chamber with filter	
	Roller			Roller	
	Latex tube injection site			Y injection site	
	Without needle			Without needle	

Blood transfusion set



1. protective cap of closure-piercing device
2. closure-piercing device
3. integral air-inlet with air filter and closure
4. fluid channel
5. drip tube
6. drip chamber
7. fluid filter
8. tubing
9. flow regulator
10. Injection site(latex tube)
11. Needle
12. protective cap of needle
13. Injection site (Y site)

Product	Components	Ref.code	Product	Components	Ref.code
	Spike with air vent	N0101		Spike with air vent	N0201
	Drip chamber with filter			Drip chamber with filter	
	Roller			Roller	
	Latex tube injection site			Y injection site	
	With needle			With needle	
	Spike with air vent	N0500		Spike with air vent	N0600
	Drip chamber with filter			Drip chamber with filter	
	Roller			Roller	
	Latex tube injection site			Y injection site	
	Without needle			Without needle	





Hitec Medical Co.,Ltd [Shanghai]

Add: No. 501, Building 2 Lane 1505,
Lian Hang Road Minhang District, 201112 Shanghai
PEOPLE'S REPUBLIC OF CHINA

Jiangsu Hitec Medical Co.,Ltd

Add:North Jiangxin road Jianghai Linkage Development zone,
Tongzhou Bay 226400 Nantong City, Jiangsu Province

Email:sheng@hitecmed.com

Mobile:00 86 134 7253 3440(Mr. Jackie Sheng)

Tel:00 86 21 64296678

www.hitecare.com