



Verifine® Safety Lancets-M Type

Less Packaging Volume, the Same Promisemed Quality and Reliability



Promisemed Verifine® Safety Lancets-- M Type

is about 30% packaging volume less than our T type to save freight fee and storage space without sacrificing premium features.



FEATURES AND BENEFITS

Safe

- The needle sterilized by irradiation is fully shielded before use
- Split-second needle retraction technology prevents used needle from being exposed
- Special design eliminates reuse

Comfortable

- Special designed tri-bevel needle tip to reduce contact surface with skin
- Double spring achieves triggering and retraction within 0.03 second
- Ergonomics design for easy hold during the puncture

Reliable

- Fully automatic assembly and testing to ensure stable quality
- Puncture depth is independent of triggering force
- Premium material to protect the needle tip when the sheath is pulled out

Simple to Use

- One-step activation
- No pre-loading required

Cost-effective

- About 30% packaging volume less than our T type to save freight fee and storage space

Environment Friendly

- Less material is used to protect environment




ORDERING INFORMATION

REF	Inside Code	Needle Gauge	Needle Length	Qty/Box	Qty/Case
SL-30-18M	SL-30-18M2	30G	1.8 mm	100	2000
SL-28-18M	SL-28-18M2	28G	1.8 mm	100	2000
SL-23-18M	SL-23-18M2	23G	1.8 mm	100	2000
SL-21-18M	SL-21-18M2	21G	1.8 mm	100	2000
SL-21-24M	SL-21-24M2	21G	2.4 mm	100	2000

※ If you have further questions or need other specifications, please contact us.

INSTRUCTIONS FOR USE



PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Cangxing Street, Cangqian Community,
Yuhang District, Hangzhou City, 311121 Zhejiang, China
info@promisemed.ca www.promisemed-global.com