

Verifine®

Safety Lancet



Your
ideal choice
for Safety
Lancet!

- **Safety**

Needle is fully shielded before and after use
Sterilized by gamma irradiation

- **Comfortable**

High quality tri-bevel with high-speed penetration
Tri-levelground cannula keeps good blood flow
The optimum ergonomics design

- **Convenient**

Pre-loaded and pressure activated device
Less pressure powder




 **Promisemed®**

Ordering Information

Catalog	Needle Size	Needle Length	Package pcs/box
SL-18-18	18G	1.8 mm	100
SL-21-12	21G	1.2 mm	100
SL-21-18	21G	1.8 mm	100
SL-21-22	21G	2.2 mm	100
SL-21-24	21G	2.4 mm	100
SL-21-28	21G	2.8 mm	100
SL-23-18	23G	1.8 mm	100
SL-23-20	23G	2.0 mm	100
SL-23-24	23G	2.4 mm	100
SL-26-18	26G	1.8 mm	100
SL-28-16	28G	1.6 mm	100
SL-28-18	28G	1.8 mm	100
SL-30-16	30G	1.6 mm	100
SL-30-18	30G	1.8 mm	100

Shelf life (sterility): 5 years old.

Test Procedure



PROMISEMED MEDICAL DEVICES INC.
 170 – 422 RICHARDS STREET,
 VANCOUVER BC V6B 2Z4, CANADA

Website: www.promisemed.ca

